

Site Number: \_\_\_\_\_

Screening ID: \_\_\_\_\_ - \_\_\_\_

First 3 Letters of First Name: \_\_\_\_\_

## A. COLLECTION INFORMATION

- Label one **10 ml** plain red top tube at **Week 2** (same tube used to collect blood for both MMF and DZB levels). At **Week 4**, label one **5 ml** plain red top tube. Do not use a barcode for this tube.
- Label two 1.8 ml Nunc cryovials with the appropriate barcode labels indicating specimen type (DZB) from an **unused** DZB barcode label sheet. Write in the first three letters of the participant's first name and the date of draw (MM/DD/YYYY) with an alcohol-proof pen. Apply the label vertically.
- Attach matching barcode labels from the **same barcode sheet** to each page of this Specimen Transmittal Form in **Section B**.
- Draw blood into the **10 ml** plain red top tube (on **Week 4** draw blood in a **5 ml** plain red top tube). Rotate tube gently and place upright in tube rack. Allow blood to clot for 15-30 minutes at room temperature (65-75° F).
- Centrifuge for 10-15 minutes.
- Split serum between two appropriately labeled 1.8 ml cryovials. Screw the tops on tightly to prevent leakage.
- Place samples upright in a 3" partitioned freezer storage box. Samples from multiple subjects may be placed in the same box.
- Freeze at -20° C.

1. Date specimen collected: \_\_\_\_\_

\_\_\_\_ / \_\_\_\_ / \_\_\_\_  
MM DD YYYY

2. For which visit, in the study sequence, is this form being completed? (*check one*)

☐ 4

Week 2

☐ 6

Week 4

☐ 99

Other

## B. SPECIMEN INFORMATION

1. Place DZB Barcode Label Here:



## C. SHIPPING INFORMATION

- Place the 3" partitioned freezer storage box into a biohazard Ziploc bag with an absorbent sheet.
- Place the yellow copy of this completed form into the outside sleeve of the bag.
- Place the bag into a large styrofoam box filled to capacity with dry ice (at least 5 lbs or 3 kg) and tape outer cardboard box securely closed.
- Affix the following **two** labels to the outside of the box: 1) Black Diamond UN1845 Dry Ice Label, 2) Diamond UN 3373 Diagnostic Specimen Label (Placed on the same side as the Black Diamond Dry Ice Label)
- Prepare and print a pre-paid airbill to FedEx all samples to **Priority Overnight**: Specimen Processing, Northwest Lipid Research Laboratories, 401 Queen Anne Avenue North, Seattle, WA 98109-4517 **Phone:** (206) 685-3327
- Ship specimens Monday – Thursday only** (except days before U.S. federal holidays).

1. Shipped By Name: \_\_\_\_\_

2. Phone #: (\_\_\_\_) \_\_\_\_ - \_\_\_\_

3. Date Shipped: \_\_\_\_ / \_\_\_\_ / \_\_\_\_  
MM DD YYYY

4. Comments: \_\_\_\_\_

## D. For TrialNet Core Lab Use Only

Sample Received?

Y

N

Date Received:

\_\_\_\_ / \_\_\_\_ / \_\_\_\_  
MM DD YYYY

Place Lab Barcode Label Here

Comments: \_\_\_\_\_

On all questions write "?" if the desired information is currently unavailable, but is being checked and will be known in future updates.  
Write "\*" if the desired information is permanently unavailable (i.e. will not be known in any future updates).