

AFRICAN AMERICAN STUDY OF KIDNEY DISEASE AND HYPERTENSION
AASK COHORT STUDY
BLOOD PRESSURE ANEROID CALIBRATION FORM # 109

The Tycos Classic Hand Aneroid device should be calibrated once every six months against a mercury sphygmomanometer. This form is completed for each Tycos Classic Hand Aneroid device at your center once every six months.

- 1. Clinical Center Number ... CC-N
2. Device number, being calibrated ... Device-num
3. Date of calibration (mm/dd/yyyy) ... / / Calibrate dt
4. Site where calibration was performed ... perform-site
5. Tech ID that is performing the calibration ... tech-ID
6. At which levels did you successfully check the calibration of the aneroid device to within +4 mmHg of the mercury sphygmomanometer (based on the average of two readings)? (0=no, 1=yes for each level.)
a. 250 mmHg ... mmHg-250
b. 160 mmHg ... mmHg-160
c. 70 mmHg ... mmHg-70

IMPORTANT: If at any of these three levels you could not check the calibration of the aneroid device to within +4 mmHg of the mercury sphygmomanometer (based on the average of two readings), you MUST return the aneroid device to Welch Allyn Medical Division following the instructions in Chapter 6 of the AASK Cohort Study Manual of Operations.

- 7. If 6 a, b, or c is answered 0=no, what date was the aneroid device mailed to Welch Allyn Medical Division? (mm/dd/yyyy) ... / / MAIL-DT
200. Date this form completed (mm/dd/yyyy) ... / / COMPL-DT
201. Certification ID of person reviewing this form ... COMPL-ID

For Clinical Center Use Only:
Certification ID of person entering this form: \_\_\_\_\_
Date Entered: \_\_\_\_/\_\_\_\_/\_\_\_\_