Archived Urine Sample Collection Form

This form is to be completed at visit 10 and 12.

ARCHIVED URINE SAMPLES FOR THE NIDDK BIOSAMPLE REPOSITORY

Collection: A freshly voided urine sample will be collected. The urine specimens will be centrifuged in 50 mL tubes at 500 g for 5 minutes as soon as possible, with volume, processing times and voiding times noted (processing times should be no longer than 20-30 minutes from the time of acquisition). Tubes will be kept in ice throughout the process. The bottom 250 μ L pellet (sometimes barely- or non-visible) will be transferred with a 1.0 mL pipette to a 1.5 mL eppendorf tube previously prepared with 750 μ L of TriReagent and inverted several times and put on ice prior to freezing at -80° C.

The remaining urine samples will then be transferred to 10 mL polypropylene (not polystyrene) Falcon culture tubes, stored in six 5 mL aliquots.

Storage Instructions: Samples are to be stored at the site (-80C) for up to four months after collection.

Shipping Instructions: Samples are to be batch-shipped quarterly on at least five pounds of dry ice. Send the freshly voided urine samples to the NIDDK Biosample Repository at Fisher Bioservices. Use pre-printed Fed Ex airbills. Do not ship on Friday.

1. LABELS:

		Check if collected	Type of Sample	Bar Code Label
		☐ fvcoll	A. Freshly Voided Urine Biosample	fvlabel
		□ pecoll	B. Freshly Voided Urine Pellet	pelabel
	Date of 0	Date of Collection:/ dtcoll		
	Voiding T	ime:	Volume:	Processing Time: proctime
Comments:				
	Shipping Instructions: Complete Shipping Manifest and pre-printed Fed Ex airbill addressed to Fisher BioServices Coorporation (NIDDK Biosample Repository). Ship samples on dry ice per guidelines provided by Fisher. CRISP Member completing this form			
Date Form Completed / / Cdidnum Date Form Completed / / Cddate Data Entry Status: Please check to indicate that the above information has been entered □				
				nformation has been entered
Primary Entered by:			: // dedate	
Secondary Entered by: Date/				