



Genetics of Kidneys in Diabetes Study

March 18, 2002

Form 222.3GK

Page 1 of 1

Genetic Specimens Mailing List

This mailing list is used whenever the GoKinD clinic ships a container of specimens to the Central Biochemistry Laboratory (CBL) for storage and subsequent genetic analysis. If a proband/relative is in End Stage Renal Disease (ESRD) indicate 1 in the column below. If the proband/parent is part of a trio indicate 2 in the column below; 1 if single; and 3 if other.

SHIPPING INSTRUCTIONS: Three (3) tubes will be collected for each relative (2 blue top vacutainers and 1 lavender top). The specimens are not to be refrigerated. Process the lavender tube as specified in chapter 8. Place all 3 tubes from a single participant into a styrofoam 6-tube mailer then put the mailer into a cardboard sleeve. Place this into a large zip-seal plastic bag before putting the mailer in the outer mailing box. Multiple mailers may be put in one shipping box. Include the mailing list and all unused labels within the outer box.

Samples **MUST ARRIVE AT CBL WITHIN 24 HOURS OF DRAW.** Deliveries are accepted Monday-Saturday.

Three copies of this form are to be distributed as follows:

- (1) **ORIGINAL COPY** Complete and place inside insulated shipping container with specimens.
Mail to: GoKinD Central Biochemistry Laboratory
 Fairview University Medical Center
 Room L275 Mayo, Memorial Bldg.
 420 Delaware Street, S.E.
 Minneapolis, MN 55455
 Telephone: (612) 273-3645
- (2) **YELLOW COPY** Send to the Coordinating Center in the weekly forms mailing.
- (3) **PINK COPY** Retain in clinic files.

Clinic Number: ___ ___ ___ **Clinic**

Specimens Shipped on: ___ / ___ / ___ **Bvshipdt**
 Month Day Year

Specimens Collected From: ___ / ___ / ___ through ___ / ___ / ___
 Month Day Year Month Day Year
Bvcolfdt **Bvcoltdt**

LYM/DNA	BARCODE	FAMILY ID NUMBER	PROBAND	PROB/REL	DATE SPECIMEN			CPT	DNA	ESRD	TRIO
			RELATIVE	INITIALS	COLLECTED			NUMBER	NUMBER	0 OR 1	CODE
			CODE	F M L	MONTH	DAY	YEAR	OF TUBES	OF TUBES		1,2,OR 3
	Bvbarcod	Bvfamily	Bvrelcode	Bvrelinit	Bvcoldt			Bvcpt	Bvdna	Bvesrd	Bvtrio
	Bvbarcod2	Bvfamily2	Bvrelcod2	Bvrelinit2	Bvcoldt2			Bvcpt2	Bvdna2	Bvesrd2	Bvtrio2
	Bvbarcod3	Bvfamily3	Bvrelcod3	Bvrelinit3	Bvcoldt3			Bvcpt3	Bvdna3	Bvesrd3	Bvtrio3
	Bvbarcod4	Bvfamily4	Bvrelcod4	Bvrelinit4	Bvcoldt4			Bvcpt4	Bvdna4	Bvesrd4	Bvtrio4