



Attention - DO NOT enter patient data on this form if the header does not contain *preprinted* HALT PKD ID number, clinical center ID, and visit number.

Participant ID: _____ *halted* Clinical Center: _____ *clinic* Date of Collection ____/____/____
dym dvd dvy
visit: _____ Form not completed misfrm

Missing Data Codes: A-Participant Refused B-Reading Not Possible C-Institutional Error

ARCHIVED BLOOD SAMPLE COLLECTION FORM

Form # 18

This form is to be completed at Visit B1, then annually. It is to be entered within two weeks of sample collection. If Missing Data Codes are used, enter code in both fields of Collection Time (e.g., BB:BB). Samples must be shipped on the day of collection to the NIDDK Biosample Repository at Fisher BioServices.

LABELS:

Affix the "Bio-serum" and "Bio-plasma" labels to this form (accession number + suffix -S for serum, suffix -P for plasma).

Affix corresponding labels on both tubes per sample.

1. Serum Sample: Collect 2 SST tubes of blood (tiger-top, 10 ml each). Gently invert 5 times, but do not shake. Allow to clot in a vertical position for 30 minutes. Centrifuge at 1300 RCF (g) for 15 minutes (within 1-2 hours of collection). Refrigerate samples. No decanting is necessary. Ship sample on the day of collection per instructions below.

2. Plasma Sample: Collect 2 PST tubes (green/grey cap, 8 ml each). Gently invert 8-10 times, but do not shake. No clotting time is necessary. Centrifuge at 1300 RCF (g) for 10 minutes (within 1-2 hours of collection). Refrigerate samples. No decanting is necessary. Ship sample on the day of collection per instructions below.

Type of Sample	Collection Time 24hr	Accession# / Bar Code Label
A. Serum Sample Label: "Bio-serum" (suffix -S)	____ : ____ <i>ser1hour ser1min</i>	<i>seraccno</i> Place Label Here
B. Plasma Sample Label: "Bio-plasma" (suffix -P)	____ : ____ <i>pla1hour pla1min</i>	<i>plaaccno</i> Place Label Here

3. Comments: _____ *commnt*

Shipping Instructions: Complete Shipping Manifest and pre-printed Fed Ex airbill addressed to Fisher BioServices Corporation (NIDDK Biosample Repository). Ship samples on cold packs, per IATA 650 guidelines, but do not allow samples to freeze. If samples are shipped on a Friday, use an airbill marked for Saturday Delivery (addressed to Heather Higgins at the nearest Fed Ex hub).

HALT PKD staff member completing this form: _____ *cmidnum* Date: ____/____/____
Month cdm Day cdd Year cdy

Data Entry Status: Please check to indicate that the above information has been entered

Primary Entered by: _____ *deidnum* Date: ____/____/____ *dem/ ded/ dey*