


PREPARATION AND SHIPPING OF VEINS

Histology Core Facility
Dr. Charles Alpers
University of Washington
Seattle

Preparation, shipping, processing




Required materials and preparations

1. Color coded and labeled 2ml tubes with solution (supplied by Alpers Lab)
2. Cryomolds, dry ice, aluminium foil.




3. Disposable blades, forceps, gloves, permanent histopen (StatmarkPen), cutting board (or disposable Petri dishes), cryotubes.



Tissue preparation

1. Following excision of vein, divide into 4 even transverse sections using clean blades.



What if ?

Q. There is less than 5mm of tissue.

A. Keep the thickness of the slices at 1.25 minimum.

If one slice of tissue (1.25mm thick) is available, preserve it in Formalin

If two slices of tissue are obtained – preserve in Formalin and Frozen OCT

If three slices of tissue are obtained – use for Formalin, Frozen OCT, RNAlater

if four slices of tissue are obtained – collect Formalin, Frozen OCT, RNAlater, and Proteomics

1. Formalin


2. Formalin , Frozen OCT

3. Formalin, Frozen OCT, RNAlater

4. Formalin, Frozen OCT, RNAlater, Proteomics




Q. I have more than 5mm of tissue.

A. Divide the available tissue into the 4 even slices.




Tissue preparation

1. Following excision of vein, divide into 4 even transverse sections using clean blades.
2. Move the blades away from each other gently cutting through the vein. Avoid pushing down and crushing the tissue.




Tissue preparation – video 1

1. Following excision of vein, divide into 4 even transverse sections using clean blades.
2. Move the blades away from each other gently cutting through the vein. Avoid pushing down and crushing the tissue. Keep the blades as close to each other as possible.



Tissue preparation for OCT Frozen – video 2

Prepare aluminium foil custom mold using cryomold.



Tissue preparation for OCT Frozen – video 3

Precool all parts of the prep on dry ice for at least 1 minute before placing the tissue slice.




Tissue preparation for OCT Frozen – video 4

Place the tissue in the middle of the aluminum mold, transfer it into vinyl cryomold and gently pack with second piece of foil.


Keep the cryomold up at all times. Try to use forceps as often as possible to avoid heat transfer from fingers.

Write necessary identifying information on the foil to match the info on the cryomold. Leave on dry ice and transfer to -20, or -80 for storage.



Tissue preparation

2. Place one segment into each of 3 different tubes.



10%
Buffered
Formalin





Proteomic
Buffer



RNAlater
Media


Formalin fixation

1. Allow to fix for 24 to 72 hours at room temperature and then replace the formalin with 70% ethanol.
 - option 1: Pipet the formalin out and pipet 70% ethanol in. Mark the tube using a permanent marker, or a new label (“70% ethanol, [date]”). Store at room temperature (RT).
 - option 2: Transfer tissue using forceps to a new tube with 70% ethanol. Store at RT.




RNAlater media

1. Allow to sit at least overnight at 4⁰ C (tissue RNA is stable at 4⁰ C up to 1 week).
2. Move the tissue to a cryotube without solution and place at -20 or -80⁰ C for longterm storage.
3. 500ml of RNAlater costs around \$300. 2ml used per sample, cost of RNAlater per sample/tube ~ \$ 1.2



RNA tissue collection




24h at fridge
4⁰ C



Long term storage
-20⁰ C or -80⁰ C,
a cryotube without solution



Where and when to send?



Formalin
fix tissue
in Ethanol
in RT

OCT on
DRY ICE

To UW of Seattle, Histology Core

1. Green (white) tube formalin fixed tissue in 70% ethanol. Send in room temperature (RT).
2. OCT frozen tissue in cryomold. Send on dry ice. Samples are sent by the Clinical Center once a month. Frozen OCT tissue must be sent overnight on **DRY ICE** only, Monday to Thursday only, with contact information forwarded to UW at least a day before. UW will confirm the readiness to receive that shipment (email or call UW before you send tissue).



Proteomic
?

RNA later
tissue on
DRY ICE


To Fisher (NIDDK Repository)

1. RNA later tissue. Send overnight on **DRY ICE** only, Monday to Thursday only, with contact information forwarded to Fisher at least a day before . Fisher will confirm the readiness to receive that shipment (email or call Fisher first before you send tissue).
2. Proteomic tissue. To be determined and TBA.



If you need to write on tubes (especially cryotubes) use StatmarkPen, or any Histopen. Do not use regular permanent marker (Sharpie)!!!

Preparation, shipping, processing



Materials for 600 pts.

2539

| | vol., case of | will need | \$ per | total | | | |
|-------------------------|---------------|-----------|--------|-------|--|----------------------------|----------------|
| Formalin | 4 l | 1x | 40 | 40 | | vwr | |
| 70 %Ethanol | 4 l | 1x | 103 | 103 | | vwr | |
| RNAlater | 0.5 l | 2x | 300 | 600 | | vwr | |
| Proteomics | 1.2 l | | ? | ? | | | |
| OCT | | 12x | 9.33 | 95 | | tissue-tek | |
| 2-methylbutane | 0.5 l | 6x | 61 | 366 | | vwr | |
| | | | | | | | |
| | | | | | | | |
| Materials | | | | | | | |
| 2ml tubes green | 500 | 2x | 46.5 | 93 | | us scientific | |
| 2ml tubes blue | 500 | 2x | 46.5 | 93 | | us scientific | |
| 2ml tubes red | 500 | 2x | 46.5 | 93 | | us scientific | |
| cryomolds | 100 | 6x | 19.5 | 117 | | tissue-tek | |
| cassettes | 600 | 1x | 100 | 100 | | vwr | |
| histopen | 1 pack of 10 | 2x | 30 | 60 | | vwr | |
| blades | 1 pack of 100 | 12x | 20 | 240 | | vwr | |
| labels laser | 3120 | 2x | 53 | 106 | | small cryo-babies by tough | |
| forceps | 1 | 6x | 3 | 18 | | vwr #25601-008 | |
| petri dishes non-vented | 480 | 1 x | 63 | 63 | | vwr #82050-570 | |
| shipping | | | | | | | |
| boxes thermosafe | 4 | 2x | 64 | 128 | | vwr #14100-432 | 25.4x25.4x17.8 |
| ups soft | | | free | | | ups | |
| cardboard box | 1 | 64 | 3.5 | 224 | | us scientific | |

dry ice

| | | |
|-----------|------|-----|
| overnight | 5 lb | 8lb |
| ups | | |
| fedex | | |

| |
|------|
| 10lb |
| |
| |

| | |
|-------|----------|
| | envelope |
| ups | |
| fedex | |

Questions?

University of Washington
Alpers Lab
HSB-Room E506
1959 NE Pacific St
Seattle, WA 98195

Contact Tomasz Wietecha, Kelly Hudkins
(206)-543-6746
tomaszw@u.washington.edu
Cell (emergencies): (425)-681-0275
(Tomasz)