| Diabetes   |   |          |                 |                                   | MF-DZB          |         | Form MMF99EL |         |              |                          |
|--|---|----------|-----------------|-----------------------------------|-----------------|---------|--------------|---------|--------------|--------------------------|
| TrialNet   |   |          | I               | ELISPOT SPECIMEN TRANSMITTAL FORM |                 |         |              |         |              | October 01, 2005         |
|  |   |          |                 |                                   |                 |         |              |         |              | Version 2<br>Page 1 of 1 |
| Site Number:   |   |          |                 | Screening ID: First 3 Lette       |                 |         |              | Letters | of First Nam | - C                      |
| A. COLLECTION INFORMATION  |   |          |                 |                                   |                 |         |              |         |              |                          |
| • I  | Label two 10 ml green top (sodium heparin) blood collection tubes with a subject identifier.  |          |                 |                                   |                 |         |              |         |              |                          |
| • I  | Label each 10 ml tube with the appropriate barcode label indicating specimen type (ELSPT) from an <b>unused</b>                                 |          |                 |                                   |                 |         |              |         |              |                          |
|  | ELSPT barcode label sheet. Write in the first three letters of the participant's first name and the date of draw (MM/DD/YYYY) with an           |          |                 |                                   |                 |         |              |         |              |                          |
| а  | alcohol-proof pen. Apply barcode label vertically.  |          |                 |                                   |                 |         |              |         |              |                          |
| • ,  | • Attach matching barcode labels from the <b>same barcode sheet</b> to <u>each page</u> of this Specimen Transmittal Form in <b>Section B</b> . |          |                 |                                   |                 |         |              |         |              |                          |
| • I  | • Draw blood into two 10 ml green top tubes and gently invert the tube 6-8 times to mix the sample. DO NOT CENTRIFUGE.                          |          |                 |                                   |                 |         |              |         |              |                          |
| • ]  | • Keep tubes at ROOM TEMPERATURE. Ship as whole blood PRIORITY OVERNIGHT the same day as blood is drawn!  |          |                 |                                   |                 |         |              |         |              |                          |
| • 9  | • SAMPLES MUST BE SHIPPED ON THE DAY OF THE BLOOD DRAW!!  |          |                 |                                   |                 |         |              |         |              |                          |
|  |   |          |                 |                                   |                 |         |              |         |              |                          |
| 1. Date specimen collected: $\frac{-\frac{1}{MM}}{DD} = \frac{1}{YYY}$ |   |          |                 |                                   |                 |         |              |         |              |                          |
| 2. F   | 2. For which visit, in the study sequence, is this form being completed? (check one)  |          |                 |                                   |                 |         |              |         |              |                          |
|  | $\square_2$   | Baseline | □ <sub>14</sub> | Month 6                           | □ <sub>23</sub> | Month 1 | 18           | □ 99    | Other        |                          |

## **B. SPECIMEN INFORMATION**

Month 3

1. Place ELSPT Barcode Label Here:



□ 29

Month 24

## C. SHIPPING INFORMATION

- Place the two 10 ml blood collection tubes of whole blood into a styrofoam tube holder with an absorbent pad.
- Place the styrofoam tube holder into a cardboard sleeve and then into a biohazard Ziploc bag.

Month 12

- Place the yellow copy of this completed form in the outside sleeve of the bag.
- Ship sample at **ROOM TEMPERATURE** in a styrofoam shipping container. Tape outer box securely closed.
- Affix the following label to the outside of the box: "Diagnostic Specimens"
- Prepare and print a pre-paid airbill to FedEx all samples **Priority Overnight** to:

**UCHSC** 

Barbara Davis Center

ATTN: Rebecca Wagner /Gottlieb Lab 1775 North Ursula Street, Room 4201U

Aurora, Co 80010 Phone: (303) 724-6804

• Ship specimens **Monday-Thursday** (except days before a U.S. federal holiday)

| 1. Shipped By Name: | <br>2. Phone #: | ()   |
|---------------------|-----------------|--|
| 3. Date Shipped:    |                 | $\frac{1}{1}$ $\frac{1}$ |
| 4. Comments:        |                 |  |
|                     |                 |  |
|                     |                 |  |

On all questions write "?" if the desired information is currently unavailable, but is being checked and will be known in future updates.

Write "\*" if the desired information is permanently unavailable (i.e. will not be known in any future updates).