

HCYA Screen Content

Homocysteine

1. tHcy value: _____ $\mu\text{mol/L}$ (DMS ACCEPTABLE RANGE 3.0 - 30.0) HCYA1

2. Specimen Code: _____ HCYA2

- A Normal
- B Hemolyzed
- C Thawed
- D Insufficient volume
- E No Sample
- F Instrument Failure
- G Technician Error
- H Other

Specify: _____ (automatic notelog)

3. Date analyzed: ____ / ____ / _____ HCYA3
mm/dd/yyyy